




PROTECTIVE FACE MASK

**N95 / FFP2
APPROVED**



EU 149:2001 + A1:2009 – FFP2 NRD



Maschera filtrante N95/FFP2 per la protezione da particelle N. repertorio: 1936777

Description:

CE&FDA Certificated

Certified N95/ FFP2 Disposable face mask
Laboratory Tested and official certificated CE and FDA
High performance filtering media with low breathing resistance
Exceptional User Comfort and Fit
Not Reusable

Features

- * Packing: boxed (4 pcs per box)
- * Size: 16x15x0.23cm
- * Layered, Non-woven Melt blown filtering media
- * Higher Electrostatic Charged to attract and protect against smaller particles
- * Increased efficiency filtering media for the highest particle protection above standard requirements
- * Premium Quality Melt Blown Media
- * FDA, CE qualified

Qualities

- * Flat fold design adds to a better fit and higher sealing around different face shapes and nationalities
- * Premium Elastane neck strap. (Maintains shape and length)
- * Hypoallergenic inner layer. (Prevents skin irritation)
- * High Quality Nose Foam (Comfortable on nose bridge)
- * Foldable Mask for easy Go

How to use:



1. Unfold the mask to cover the nose and jaw



2. Gently press the setting line to fit the nose



3. Hang earloop on both ears



4. Adjust to ensure the edge of the mask sticks to the face



Certificate of Conformity

CERTIFICATE OF CONFORMITY

Certificate No.: CCT20031201TCS

Applicant.....: ██████████ology Co., Ltd.
██████████7, Sungang Village, Baogang Road, Sungang Street,
Luohu District, Shenzhen,518100,CN

Manufacturer.....: ██████████ology Co., Ltd.
██████████Sungang Village, Baogang Road, Sungang Street,
Luohu District, Shenzhen,518100,CN

Product.....: **N95 face mask**

Trade Name.....: **STARBUK**

Model(s): ██████████K-
██████████MK-NX04


It is only valid in connection with the test report number CCT20031201TRS Identification of regulation/standards

EN 149: 2001+A1:2009

The EUT described above has been tested by us with the listed standards and found in compliance with the council PPE Directive **EU 2016/425**. It is possible to use CE marking to demonstrate the compliance with this PPE Directive.

This certificate of conformity is based on a single evaluation of the submitted sample(s) of the above mentioned product. It does not imply an assessment of the whole production and other relevant directives have to be observed.




Toby/Senior Manager
Mar. 12, 2020



Shenzhen CCT Testing Technology Co., Ltd.
8th Floor, Area I, Building 1, Hanhaida Science and Technology Innovation Park, Yulv Community,
Guangming New District, Shenzhen, Guangdong, China
Web: www.cct-prc.com Tel: 400-8788-298 Fax: 0755-33609691 Email: cct@cct-prc.com