



PROTECTIVE FACE MASK

**KN95 / FFP2
APPROVED**



EU 149:2001 + A1:2009 – FFP2 NRD



Maschera filtrante KN95/FFP2 per la protezione da particelle - N. repertorio: 1936939

Description:

10Pack of KN95 Face Masks, Multi-Ply Cotton Filter Medical Sanitary for Dust, Germ Protection:

KN95 Face Masks, Multi-Ply Cotton Filter Medical Sanitary for Dust, Germ Protection and Personal Health Care

Protective face masks made from non-woven material and feature elastic earloops for comfort. Bacterial Filtration Efficiency (BFE) is the effectiveness of a material to filter bacteria of specified particle size, 3.0 μm (micron) $\pm \mu\text{m}$. This expressed as a percentage of a quantity that does not pass through the material.

Features

- * Packing: boxed
- * Size: 6.3"x4.72"x2.75" (16x12x7cm)
- * PFE95%, with a $\geq 95\%$ protective efficiency * Anti PM2.5
- * Multi-ply cotton filter, that is comfortable and Hypoallergenic for wearing
- * With ear loops for wearing
- * Fluid Resistance/Flammability/Needle Filter Cotton
- * FDA, CE qualified

Notes:

Masks contain an aluminium nose piece and should be stored in a cool and dry place, avoiding excessive heat (40°C / 104°F).

Masks should be disposed of after every use.

HOW TO USE

- * Pull and extend the mask with two hands and be sure that the nose clip is upward
- * Hold the mask onto face and pull elastic bands over ears
- * Press and adjust nose clip to fit and form the shape of nose bridge
- * Discard after every use



Multi-layer Protection

- Multi-layer Protection
- Up to 95% filtration efficiency
- Air pollution resistance (PM 2.5)
- Medical-grade filter protects against harmful particles
- Comfortable 3D flex design
- Reusable & eco-friendly
- Suitable for adults and children

FDA ISO CE EN149



Certificate – 證明書 – Сертификат – Certificat – 증명서 – certificado

Certificate of Conformity

Certificate No. : STA20200313008EC

Certificate's Holder : ~~XXXXXXXXXXXXXXXXXX~~ CO., LTD

Address : ~~XXXXXXXXXXXXXXXXXX~~ Jishui County, Ji'an City,
Jiangxi Province, China.

Manufacturer : ~~XXXXXXXXXXXXXXXXXX~~ CO., LTD

Address : ~~XXXXXXXXXXXXXXXXXX~~ County, Ji'an City,
Jiangxi Province, China.

Product Name : KN95 Face Mask

Model(s) : ZC9595

Trade Mark : ~~XXXXXXXXXX~~

Test Report Number(s) : SCC(20)-60068S-PPE

Test Standard : EN 149:2001+A1:2009

The submitted products have been tested by us with the above standards and found in compliance with the following European Directives:

PPE Directive : 2016/425/EU

Note:

This certificate is only valid for the equipment and configuration described, and in conjunction with the test data detailed above. The CE mark as show joined can be used, under the responsibility of the manufacturer or the importer, after completion of the CE Declaration of Conformity and in accordance with the above directives.

Issuing Date: March 18, 2020



Shenzhen STA Testing Co., Ltd.
Room 204, Fuyong chamber of commerce information building, Road 6348, Baoan District, Shenzhen, China
☎ : www.xunda-test.com ☎ : 86-755-29985395 ✉ : sta-test@xunda-test.com